General Immunofluorescence Assay Method

PRODUCT INSERT

INTENDED USE

The components listed below may be used for detection of antibodies in human serum by indirect immunofluorescence.

SUMMARY AND EXPLANATION

These instructions for use provide a general protocol for use of substrate slides and conjugates appropriate for indirect immunofluorescence studies. Autoantibodies may be visualized on a variety of tissue section and cell culture substrates using a FITC conjugate system. The principle of the procedure is described in the following section. Individual laboratories may use this product insert as a guide but should perform internal validation of internal protocols using these materials.

PRINCIPLES OF PROCEDURE

The components listed below may be utilized in indirect immunofluorescence protocols. In such a protocol, patient sera is incubated on a substrate slide to allow binding of antibodies. Any antibodies not bound are removed by rinsing. Bound antibodies of the IgA, IgG and/or IgM class, as determined by the conjugate used, are detected by incubation of the substrate with fluorescein-labeled, anti-human conjugate. Reactions are observed under a fluorescence microscope equipped with appropriate filters. The presence of antibodies may be demonstrated by an apple green fluorescence of specific histological structures in the substrate. The titer may be determined by serial dilution of the specimen. The reciprocal of the highest dilution giving a positive reaction is the titer.¹

PRODUCT INFORMATION

Storage and preparation

Store all reagents at 2-8°C. Reagents are ready for use after equilibration to room temperature.

Materials provided

These instructions are provided for use with individual components, including the product codes listed below. Substrate slides are listed alongside the conjugate appropriate for detection of autoantibodies on the particular substrate.

Substrate Slide		Conjugate		Positive Control		Specimen Dilution	Sample Incubation Time*
[REF]	Description	[REF]	Description	[REF]	Description		
2123	Primate Adrenal Gland Slide, 6 well	2099	Anti-human IgG FITC with Evans Blue Counterstain, 5ml (primate-adsorbed)	2279	Adrenal Antibody Positive Control	1:4	30 minutes
2124	Primate Salivary Gland Slide, 6 well	2099	Anti-human IgG FITC with Evans Blue Counterstain, 5ml (primate-adsorbed)	2201	ANA Homogeneous Control	1:10	30 minutes
2125-4	Primate Ovary Slide, 4 well	2099	Anti-human IgG FITC with Evans Blue Counterstain, 5ml (primate-adsorbed)	2206	Steroidal Cell Antibody Positive Control	1:10	3 hours
2127	Primate pituitary slide, 6 well	2099 or 2100	Anti-human IgG FITC with Evans Blue Counterstain, 5ml (primate-adsorbed)	2201	ANA Homogeneous control	1:10	30 minutes
2128	Primate Cerebellum Slide, 6 well	2099	Anti-human IgG FITC with Evans Blue Counterstain, 5ml (primate-adsorbed)	2280	Hu Antibody Positive Control	1:10	30 minutes
2134	Primate Sciatic Nerve Slide, 6 well	2140	Goat antihuman IgM FITC with Evans Blue Counterstain, 5ml	2209	MAG Antibody Positive Control	1:10	30 minutes

EN

2147	Primate Split Skin Slide, 6 well	2099	Anti-human IgG FITC with Evans Blue Counterstain, 5ml (primate-adsorbed)	2213 and/or 2217	Intercellular Antibody Positive Control and/or BMZ Antibody Positive Control	1:10	30 minutes
2148	Rat Kidney/Stomach Slide, 8 well	2100	Goat anti-human IgG FITC with Evans Blue Counterstain, 5ml (primate-adsorbed)	2201, 2210	ANA (Homo) Positive Control ,AMA Positive Control	1:10	30 minutes
2149	Primate Nerve/Bovine Nerve Slide, 6 well	2140	Goat antihuman IgM FITC with Evans Blue Counterstain, 5ml	2209	MAG Antibody Positive Control	1:10	30 minutes
2155	Primate Esophagus Slide, 6 well	2100 (or 2099)	Goat anti-human IgG FITC with Evans Blue Counterstain, 5ml (primate-adsorbed)	2213 and/or 2217	Intercellular Antibody Positive Control and/or BMZ Antibody Positive Control	1:10	30 minutes
2156	Transitional Epithelium slide 6 well	2100	Goat anti-human IgG FITC with Evans blue counterstain	2241	PNP positive Control	1:10	30 minutes
2157	Primate Heart Muscle Slide, 6 well	2099	Anti-human IgG FITC with Evans Blue Counterstain, 5ml (primate-adsorbed)	2235	Heart/ Skeletal Muscle Antibody Positive Control	1:10	30 minutes
2158	Primate Skeletal Muscle Slide, 6 well	2099	Anti-human IgG FITC with Evans Blue Counterstain, 5ml (primate-adsorbed)	2235	Heart/ Skeletal Muscle Antibody Positive Control	1:10	30 minutes
2169	Mouse Stomach Slide, 6 well	2100	Goat anti-human IgG FITC with Evans blue counterstain	2201	ANA Homogeneous control	1:10	30 minutes
2172	Rat Skeletal muscle, 6 well	2100 (or 2099)	Goat anti-human IgG FITC with Evans Blue Counterstain, 5ml (primate-adsorbed)	2201	ANA Homogeneous control	1:10	30 minutes
2180	Primate Thyroid Slide, 6 well	2099	Anti-human IgG FITC with Evans Blue Counterstain, 5ml (primate-adsorbed)	2239	Thyroid/ Microsomal Antibody Positive Control	1:10	30 minutes

* The second incubation (conjugate incubation) in all cases is 30 minutes.

Optional components

[REF]	Description	U/M
2099	Primate-adsorbed goat anti-human IgG FITC with Evans blue counterstain	5ml
2100	Goat anti-human IgG FITC with Evans blue counterstain	5ml
2100x	Goat anti-human IgG FITC without Evans blue counterstain	5ml
2100x-8	Goat anti-human IgG FITC	8ml
2100-15	Goat anti-human IgG FITC with Evans blue counterstain	15ml
2100-15x	Goat anti-human IgG FITC without Evans blue counterstain	15ml
2107	Goat anti-human IgA FITC with Evans blue counterstain	5ml
2107x	Goat anti-human IgA FITC without Evans blue counterstain	5ml
2113	Goat anti-human polyvalent FITC (IgG/IgA) with Evans blue counterstain	5ml
2113x	Goat anti-human polyvalent FITC (IgG/IgA) without Evans blue counterstain	5ml
2113-15	Goat anti-human polyvalent FITC (IgG/IgA) with Evans blue counterstain	15ml
2130	Goat anti-human polyvalent FITC (IgG/IgA/IgM) with Evans blue counterstain	5ml
2130x	Goat anti-human polyvalent FITC (IgG/IgA/IgM) without Evans blue counterstain	5ml
2130-15	Goat anti-human polyvalent FITC (IgG/IgA/IgM) with Evans blue counterstain	15ml

EN

2130-15x	Goat anti-human polyvalent FITC (IgG/IgA/IgM) without Evans blue counterstain	15ml
2140	Goat anti-human IgM FITC with Evans blue counterstain	5ml
2140x	Goat anti-human IgM FITC without Evans blue counterstain	5ml
2200	IFA Negative Control	0.5ml
2201	ANA Homogeneous Positive Control	0.5ml
2206	anti-Steroidal Cell Positive Control	0.5ml
2209	anti-MAG Positive Control	0.5ml
2210	AMA Positive Control	0.5ml
2210-1	Low titer AMA Positive Control	0.5ml
2210-50	AMA Positive Control	50ml
2211	ASMA Positive Control	0.5ml
2212	AGPA Positive Control	0.5ml
2213	anti-IC Positive Control	0.5ml
2214	anti-IC Positive Control (<i>pemphigus vulgaris</i>)	0.5ml
2216	anti-IC Positive Control (<i>pemphigus foliaceus</i>)	0.5ml
2217	anti-BMZ Positive Control (pemphigoid)	0.5ml
2235	anti-Heart/Skeletal Muscle Positive Control	0.5ml
2239	anti-Thyroid/Microsomal Positive Control	0.5ml
2243	IgG FITC Conjugate with Evans Blue	60 ml
2261	Ribosomal Positive Control	0.5ml
2279	anti-Adrenal Positive Control	0.5ml
2280	anti-Hu Positive Control	0.5ml
2281	anti-Yo (PCA-1) Positive Control	0.5ml
2301	PBS	for 1 liter
2302	Sample Diluent	60ml
2500	Microscope Slide Coverglass (24x60mm)	box of 12
2505	Mounting Medium Dropper Vial	5ml
2506	Mounting Medium Dropper Vial	60ml
2510*	Evans Blue Counterstain	1.0ml

Symbols used on labels:

- [LOT] Lot number
- [REF] Catalog number
- [IVD] In vitro diagnostic use
- e Use by
- t Storage temperature
- l Consult instructions for use
- s Number of tests
- M Manufacturer
-  Date of Manufacture

* Danger. May cause cancer. Obtain special instructions before use. Do not handle until all safety precautions have been read and understood. IF exposed or concerned: Get medical advice/attention. Wear protective gloves, protective clothing, eye protection, face protection. Store locked up. Dispose of contents/container to comply with local, state and federal regulations.

Material required but not provided

- Fluorescence microscope
- Micropipette or Pasteur pipette
- Serological pipettes

EN

- Staining dish (e.g. Coplin jar)
- Small test tubes (e.g. 13 x 75 mm) and test tube rack
- Distilled or deionized water
- 1 liter container
- Wash bottle
- Paper towels
- Incubation chamber

WARNINGS AND PRECAUTIONS

For *in vitro* Diagnostic Use. All human derived components used have been tested for HbsAg, HCV, HIV-1 and 2 and HTLV-I and found negative by FDA required tests. All human serum specimens and human derived products should be treated as potentially hazardous, regardless of their origin. Follow good laboratory practices in storing, dispensing and disposing of these materials.²

WARNING - Sodium azide (NaN₃) may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal of liquids, flush with large volumes of water to prevent azide buildup. Sodium azide may be toxic if ingested. If ingested, report incident immediately to laboratory director or poison control center.

SPECIMEN COLLECTION AND PREPARATION

Only serum specimens should be used for this procedure. Grossly hemolyzed, lipemic or microbially contaminated specimens may interfere with the performance of this test and should not be used. Store specimens at 2-8° C for no longer than one week. For longer storage, serum should be frozen at -20° C. Avoid repeated freezing and thawing of samples.

PROCEDURE

Test Method

A. Screening

1. Dilute each patient serum with the Buffered Diluent according to the dilution recommended in the Materials Provided section. 1:4 (0.1 ml serum + 0.3 ml diluent) or 1:10 (10 µl serum + 90 µl diluent). Do not dilute Positive or Negative Controls to be used. Save the undiluted sera to determine antibody titers if screening tests are found to be positive.
2. Allow pouches containing substrate slides to equilibrate to room temperature for 10-15 minutes. Carefully remove the slides without touching the substrate.
3. Label the slides and place them in an incubation chamber lined with paper towels moistened with water to prevent drying.
4. Invert dropper vial and gently squeeze to apply 1 drop (or enough to assure well coverage) of the Negative Control to well #1. If applicable, apply 1 drop of Positive Control to well #2. Avoid overfilling the wells.
5. Using a micropipette or Pasteur Pipette, apply 1 drop of patient's diluted serum (or enough to assure well coverage) to the other wells. Avoid overfilling the wells.
6. Place the lid on the incubation chamber and incubate slides at room temperature according to the Sample Incubation Time recommended in the Materials Provided section.
7. Remove a slide from the incubation chamber. Hold slide at tab end and rinse gently with approximately 10 ml PBS using a pipette, or rinse slide in beaker filled with PBS. Do not use wash bottle. Transfer slide immediately into Coplin jar and wash 10 minutes. Repeat process with all remaining slides.
8. Remove slide(s) from Coplin jar. Blot the edge of the slide on a paper towel to remove excess PBS. Place the slide in the incubation chamber. Immediately invert the Conjugate dropper vial and gently squeeze to apply 1 drop to each well.
9. Repeat steps **7 and 8** for each slide.
10. Replace the lid on the incubation chamber. Incubate 30 minutes at room temperature.
11. Remove a slide from incubator. Hold the slide at the tab end and dip the slide in a beaker containing PBS to remove excess conjugate. Place slide(s) in a staining dish filled with

EN

PBS for 10 minutes. If optional conjugate without counterstain is used (see optional components in Materials Provided Section), 2-3 drops of Evans Blue counterstain may be added to the final wash. Repeat for the remaining slides. NOTE: Improper washing may lead to increased background fluorescence.

12. Remove a slide from the staining dish. Blot the edge of the slide on a paper towel to remove excess PBS. **To prevent slide from drying, proceed immediately with next step while slide is still wet.**
13. Mount the coverslip by applying **3 drops** of Mounting Medium evenly on the coverslip and place coverslip over slide. Avoid applying undue pressure and prevent lateral movement of the coverslip.
14. Repeat steps 12 and 13 for each slide.
15. Examine for specific fluorescence under a fluorescence microscope at a magnification of 200x or greater.

Slides may be read as soon as prepared. However, because of the presence of antifading agent in the mounting medium, no significant loss of staining intensity occurs if reading is delayed for up to 48 hours. Slides should be stored in the dark at 2-8° C.

B. Endpoint Determination (titration)

A serum positive in the screening test may be further tested following steps 5 through 13 to determine the titer. Make serial two-fold dilutions starting at the initial screening dilution. The reciprocal of the highest dilution producing a positive reaction is the titer. Tables for standard screening dilutions are provided below:

Preparation of Serial Dilutions Starting at 1:4

Number six tubes 1 through 6. Add 0.3 ml of Sample Diluent to tube 1 and 0.2 ml to tubes 2 through 6. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

Tubes	1	2	3	4	5	6
Serum	0.1 ml					
	+					
Buffered Diluent	0.3 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
		↗	↗	↗	↗	↗
Transfer		0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
Final dilution	1:4	1:8	1:16	1:32	1:64	1:128 etc.

Preparation of Serial Dilutions Starting at 1:10

Number six tubes 1 through 6. Add 0.9 ml of Sample Diluent to tube 1 and 0.2 ml to tubes 2 through 6. Pipette 0.1 ml of undiluted serum to tube 1 and mix thoroughly. Transfer 0.2 ml from tube 1 to tube 2 and mix thoroughly. Continue transferring 0.2 ml from one tube to the next after mixing to yield the dilutions depicted in the following table:

Tubes	1	2	3	4	5	6
Serum	0.1 ml					
	+					
Buffered Diluent	0.9 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
		↗	↗	↗	↗	↗
Transfer		0.2 ml	0.2 ml	0.2 ml	0.2 ml	0.2 ml
Final dilution	1:10	1:20	1:40	1:80	1:160	1:320 etc.

QUALITY CONTROL

EN

Both a Positive and Negative Control should be included with each test run. The Negative Control should show no specific fluorescence. The laboratory must determine the appropriate performance characteristics of the Positive Control used.

If expected results are not obtained, the run should be repeated. If inadequate results continue to occur with the controls, these may be due to:

- Turbidity. Discard and use another control
- Problems with the optical system of the fluorescence microscope. These may include improper alignment, bulb beyond useful life expectancy, etc.
- Allowing the slide to dry during the procedure.

INTERPRETATION OF RESULTS

Immunofluorescence assay results are generally reported as negative with a titer less than the screening dilution, positive with a titer greater than or equal to the screening dilution or, preferably, positive with specific endpoint titer. Laboratories must validate protocols for use of these materials and train technicians to read specific fluorescence patterns appropriate to the materials used and testing performed.

LIMITATION OF THE PROCEDURE

Results of all positive indirect immunofluorescence tests along with the results of other laboratory tests and the clinical condition of the patient when making a diagnosis.

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REFERENCES

1. Beutner EH, Kumar V, Krasny SA and Chorzelski TP. Defined immunofluorescence in immunodermatology. In "Immunopathology of the Skin", Beutner EH, Chorzelski TP and Kumar V, Eds, John Wiley and Sons, New York, 3rd Ed, 3-40, 1987.
2. Biosafety in Microbiological and Biomedical Laboratories. Centers for Disease Control, National Institutes of Health, 1993 [HHS Pub. No. (CDC) 93-8395].

For technical assistance please contact:

 **IMMCO Diagnostics, Inc.**
10 Earhart Drive
Buffalo, NY 14221-7078 USA
Telephone: (716) 691-0091
Fax: (716) 691-0466
Toll Free USA/Canada: 1-800-537-TEST
E-Mail: info@immco.com

or your local product distributor

EMERGO Europe.
Westervoortsedijk 60
6827 AT Arnhem
The Netherlands