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REFERENCER

- Lieberman J: *The specificity and nature of serum-angiotensin-converting enzyme (serum ACE) elevations in sarcoidosis*. Ann NY Acad Sci 278:488, 1976.
- Sleffkin AD, Parsons GH, Patwell SW, Hollinger MA: *The value of serial serum angiotensin converting enzyme determinations in hospitalized patients with lung disease*. Am J Med Sci 288:200, 1984.
- Alhenc-Golas F, Weare JA, Johnson RL Jr, Erdos EG: *Measurement of human converting enzyme level by direct radioimmunoassay*. J Lab Clin Med 101:83, 1983.
- Lanzillo JJ, Fanburg SL: *Development of competitive enzyme immunoassays for human serum angiotensin-1 converting enzyme: A comparison of four assay configurations*. Anal Biochem 126:156, 1982.
- Holmquist B, Bunning P, Riordan JF: *A continuous spectrophotometric assay for angiotensin-converting enzyme*. Anal biochem 95:540, 1979.
- Maguire GA, Price CP: *A continuous monitoring spectrophotometric method for the measurement of angiotensin-converting enzyme in human serum*. Ann Clin Biochem 22:204, 1985.
- Harjanne A: *Automated kinetic determination of angiotensin-converting enzyme in serum*. Clin Chem 30:901, 1984.
- Waeber B, Brunner HP, Brunner DB, Curtet AL, Turini GA, Gavras H: *Discrepancy between antihypertensive effect and angiotensin converting enzyme inhibition by captopril*. Hypertension 2:236, 1980.
- Edwards CRW, Padfield PL: *Angiotensin-converting enzyme inhibitors: Past, present and bright future*. The Lancet 1:30, 1985.
- Lieberman J: *Elevation of serum angiotensin-converting enzyme (ACE) level in sarcoidosis*. Am J Med 59:365, 1975.
- Corot C, Hentsch AM, Curtelin L: *Interaction of Gadolinium Complexes with Metal Dependent Biological Systems*. Investigative Radiology 29:2 164-167, 1994.
- Corot C, Idee, JM, Hentsch AM, Santus R, Mallet C, Goulas V, Bonnemain B, Meyer D : *Structure-Activity Relationship of Macrocyclic and Linear Gadolinium Chelates: Investigation of Transmetalation effect on the Zinc dependent Metallopeptidase Angiotensin-Converting Enzyme*. JMIRI 8:3, 696-702, 1998.
- Proctor KAS, Rao LV, Roberts WL: *Gadolinium Magnetic Resonance Contrast Agents Produce Analytic Interference in Multiple Serum Assays*. Am J Clin Pathol 121 282-292 2004.
- Jalil JE, Ocaranza MP, Pidodo AM, Jalil R: *Reproducibility of plasma angiotensin-converting enzyme activity in human subjects determined by fluorimetry with Z-phenylalanine-histidyl-leucine as substrate*. J Lab Clin Med. 1999;133:501-506.
- Domingo R, Sturrock ED, Collins M; *ACE Activity and Endurance Performance during the South African Ironman Triathlons*. Int J Sports Med. 2013;34:402-408.

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ACE-reagens

eIFU-indikator

Trinity Biotech plc,
IDA Business Park
Southern Cross Road
Bray, Co. Wicklow
Ireland
Tel. 353 1 2769800
Fax. 353 1 2769888
www.trinitybiotech.com

USA ENQUIRIES:
Trinity Biotech USA
2823 Girls Road,
Jamestown, NY 14701.
Tel. 1800 325 3424
Fax. 716-488-1990