DOA URINE TOXICOLOGY LIQUID CONTROLS

GC/MS (Low Opiate), GC/MS-H (High Opiate),
GC/MS Mix,
CONFIRM 4 (SAMHSA guidelines),
CONFIRM 5 (SAMHSA guidelines Effective Date: February 1, 2024),
PAIN MGMT (High and Low),
IMMUNOASSAY URINE TOXICOLOGY (Low, 2X, Negative, Positive, Fentanyl)

The DOA Urine Toxicology Liquid Controls are prepared from human based urine, available as a negative and at various constituent target levels to monitor the performance of qualitative or quantitative procedures for the detection of drugs in urine.

INTENDED USE

The DOA Urine Toxicology Liquid Controls are In-Vitro Diagnostic (IVD) devices, for laboratory use only, which are intended for use as quality control urine to monitor the precision of laboratory urine toxicology testing procedures for the analytes listed in the package insert.

SUMMARY AND EXPLANATION

The DEA exempt DOA Urine Toxicology Liquid Controls product line is manufactured using a human based matrix that has been stabilized to ensure that the product will be viable until the date of expiration. Positive controls are spiked with reference drug standards and/or appropriate metabolites that have been obtained from ISO certified manufacturers. Standards are certified by the manufacturers to be at least 98% minimum purity. Specific gravity, pH, and creatinine fall within the limits of normal human urine.

DESCRIPTION

Each bottle contains stabilized human based urine. Positive control urines have been spiked with authentic reference drug standards and/or appropriate metabolites. Negative control urines are certified negative by a combination of immunoassay, GC/MS and/or LC/MS for the constituents listed on our target sheets. They should be treated as any "unknown" specimen while following the specific protocol of the assay being used. *This product is intended to be used by health care professionals as an integral part of good laboratory practices.*

STORAGE & STABILITY - Please refer to Technical Note for detailed instructions.

Unopened:

- A. The controls are stable until the expiration date when stored at -10° to -20° C and protected from light.
- B. The controls are stable until the expiration date when stored at 2°-8°C, however Oxazepam and/or 6-MAM is stable for only 6 months.

After Opening:

- A. The controls are stable for 6 months or until the expiration date, whichever comes first, when stored at -10° to -20° C.
(Controls can be aliquoted and frozen)
- B. The controls are stable for 31 days or until the expiration date, whichever comes first, when stored tightly capped at 2°-8° C.
- C. Thaw controls as needed; allow to come to room temperature followed by gentle swirling before use.

Note: If end-user desires, the control can safely undergo 5 thaw/freeze cycles without affecting performance. The control can be thawed and returned to freezer to extend shelf life to 6 months. Thaw controls as needed by allowing to come to room temperature followed by gentle swirling before use.

PROCEDURE

- A. Allow controls to come to room temperature followed by gentle swirling or inversion before use. **DO NOT SHAKE.**
- B. Pipette an appropriate aliquot of DOA Urine Toxicology Liquid Control as required by the confirmation method.

Expected Results

The positive DOA Urine Toxicology Liquid Controls must test positive. The negative control must test negative. Nova Century Scientific, Inc. will (upon request), supply assay values derived from our contract assay laboratories and customer base on a particular lot of control material.

PRECAUTIONS

For In Vitro Diagnostic Use Only. Please read the entire package insert before using the DOA Urine Toxicology Liquid Controls. Please use the same safety precautions you would use for processing any "unknown" urine sample containing potentially infectious biological material. Protect product from exposure to direct sunlight. Contains sodium azide: To prevent formation of explosive metal azides dispose of waste by flushing with copious amounts of water or according to local governing regulations.

Do not use beyond the expiration date.

LIMITATIONS OF PROCEDURE

- The GC/MS, confirm 4 and 5 and Pain MGMT lines of controls are meant to be used to validate the performance of GC/MS confirmation methods.
- The Immunoassay Urine Toxicology Controls are meant to be used to validate the performance of immunoassay drug screening methods.
- Consult test manufacturer's instructions when using this product; changes in reagents, sample requirement, or methodology may effect test results.
- Although target values are provided with the DOA Urine Toxicology Liquid Controls, each laboratory should run these controls as unknowns in order to establish "in-house" assay values for them.
- *This product is not meant to be used as a standard or calibrator.*

DOA URINE TOXICOLOGY LIQUID CONTROLS, OXAZEPAM STABILITY: Oxazepam has known stability problems in urine stored refrigerated, our studies indicates that Oxazepam will deteriorate when stored refrigerated for longer than 6 months.

DOA URINE TOXICOLOGY LIQUID CONTROLS, THC STABILITY

DOA Urine Toxicology Liquid Controls are stable for the length of time under the storage conditions stated in the package insert. In spite of this fact, under certain conditions, there may be observed a gradual decline in THC levels, over time, from continuous use of a single bottle of control material. This drop in THC values may occur from any THC sample (i.e. calibrators, controls, and samples).

The apparent loss of THC most often occurs from handling and not from product instability. It is well known that THC-COOH binds to surfaces, especially certain plastics. In order to minimize this adsorption loss, we recommend the following when handling any sample (including DOA Urine Toxicology Liquid Controls) which may contain THC:

1. It is preferable to use glass pipettes or pour controls into sample cups. As an alternate, pipettors with disposable plastic tips may be used. Soft plastic transfer pipettes should be avoided.
2. Do not rinse the pipette back and forth into the sample.
3. Sample volume to surface area ratio should be as high as possible (i.e. when transferring, sample containers should be filled as much as possible with sample). Avoid rough surface plastic containers.
4. When pipetting, immerse the pipette tip as little as possible into the sample solution.
5. Do not return any unused material back into the original sample. These same guidelines should also be followed when aliquoting a control (or sample) for future use.

REFERENCES:

1. Blanc JA, Manneh VA, et al. Adsorption losses from urine-based cannabinoid calibrators during routine use. Clin Chem 1993; 39:1705-1712
2. Roth KDW, Siegel NA, et al. Investigation of the effects of solution composition and container material type on the loss of 11-nor-delta 9-THC-9-carboxylic acid. J Anal Tox 1996; 20:291-300

DOA Urine Toxicology Liquid Control

GC/MS Line Target Values (ng/mL)

CONSTITUENTS	Cutoff -25%	Cutoff	Cutoff +25%
SAMHSA MANDATED			
Delta-9-THC-COOH	11	15	19
Benzoylecgonine	113	150	188
Phencyclidine (PCP)	19	25	31
Codeine (Low Opiate)	225	300	375
Codeine (High Opiate)	1500	2000	2500
*Total Morphine (Low	225	300	375
*Total Morphine (High	1500	2000	2500
**6-Monoacetylmorphine	7.5	10	12.5
d-Amphetamine	375	500	625
d-Methamphetamine	375	500	625
NON-MANDATED			
Secobarbital	225	300	375
Phenobarbital	225	300	375
Oxazepam	225	300	375
Methadone	225	300	375
Methaqualone	225	300	375
Propoxyphene	225	300	375
Butalbital	225	300	375

* Contains Morphine-3-Glucuronide and Free Morphine (2:1)

**** 6-MAM Only in High Opiate (-H) Control**

Confirm 4 Line Target Values (ng/mL)

CONSTITUENTS	Cutoff -60%	Cutoff-25%	Cutoff	Cutoff +25%
SAMHSA MANDATED				
Delta-9-THC-COOH	6	11	15	18.75
Benzoylecgonine	40	75	100	125
Codeine	800	1500	2000	2500
*Total Morphine	800	1500	2000	2500
6-Monoacetylmorphine	4	7.5	10	12.5
Hydrocodone	40	75	100	125
Hydromorphone	40	75	100	125
Oxycodone	40	75	100	125
Oxymorphone	40	75	100	125
Phencyclidine (PCP)	10	19	25	31.25
d-Methamphetamine	100	187.5	250	312.5
d-Amphetamine	100	187.5	250	312.5
MDMA	100	187.5	250	312.5
MDA	100	187.5	250	312.5
MDEA	100	187.5	250	312.5

* Contains Morphine-3-Glucuronide and Free Morphine (2:1)

Confirm 5 Line Target Values (ng/mL)

CONSTITUENTS	Cutoff -60%	Cutoff	Cutoff +25%
Delta-9-THC-COOH	6	15	18.75
Benzoylecgonine	40	100	125
Codeine	800	2000	2500
*Total Morphine	1600	4000	5000
6-Monoacetylmorphine	4	10	12.5
Hydrocodone	40	100	125
Hydromorphone	40	100	125
Oxycodone	40	100	125
Oxymorphone	40	100	125
Phencyclidine (PCP)	10	25	31.25
d-Methamphetamine	100	250	312.5
d-Amphetamine	100	250	312.5
MDMA	100	250	312.5
MDA	100	250	312.5
MDEA	100	250	312.5

* Contains Morphine-3-Glucuronide and Free Morphine (2:1)

Pain MGMT Line Target Values (ng/mL)

CONSTITUENTS	18005116 Low	18005117 High
Codeine	120	375
Morphine	120	375
Hydrocodone	120	375
Hydromorphone	120	375
Oxycodone	120	375
Oxymorphone	120	375
Methadone	40	125
EDDP	40	125
Fentanyl	0.8	2.5
Nor fentanyl	0.8	2.5
Cotinine	80	250
Buprenorphine	4	12.5
Norbuprenorphine	4	12.5
Naloxone	10	31.25

Immunoassay Urine Toxicology Target Values (ng/mL)

CONSTITUENTS	050MO	2XMO	2XMOB	FENTANY
SAMHSA MANDATED				
Delta-9-THC-COOH	75	100	100	--
Benzoyllecgonine	450	600	600	--
Phencyclidine (PCP)	37.5	50	50	--
Morphine (Low Opiate)	450	600	600	--
d-Amphetamine	1500	2000	2000	--
d-Methamphetamine	1500	2000	2000	--
NON-MANDATED				
Secobarbital	450	600	600	--
Oxazepam	450	600	600	--
Methadone	450	600	600	--
Methaqualone	450	600	600	--
Propoxyphene	450	600	600	--
Nortriptyline	1500	2000	2000	--
MDMA*	750	1000	1000	--
Oxycodone	150	200	200	--
Buprenorphine	--	--	20	--
Norfentanyl	--	--	--	50

* 3,4 Methylenedioxyethyl-amphetamine

GC/MS Mix Target Values (ng/mL)

CONSTITUENTS	18002988
THC Glucuronide*	151
Benzoyllecgonine	1000
d-Methamphetamine	1500
d-Amphetamine	1500
MDMA	1000
MDA	1000
MDEA	1000

ORDERING INFORMATION:

CATALOG #	DESCRIPTION	SIZE	CATALOG #	DESCRIPTION	SIZE
GC/MS LOW OPIATE			Pain MGMT		
DOAC-19332001	GC/MS Cutoff -25%	50 mL	DOAC-18005116	Low Opiate	20 mL
DOAC-19330000	GC/MS Cutoff	50 mL	DOAC-18005117	High Opiate	20 mL
DOAC-19332501	GC/MS Cutoff +25%	50 mL			
GC/MS HIGH OPIATE			Immunoassay Urine Toxicology		
DOAC-19572251	GC/MS-H Cutoff -25%	50 mL	DOAC-050MO	Low Positive, Cutoff +50% with MDMA and Oxycodone	5 mL
DOAC-19572502	GC/MS-H Cutoff	50 mL	DOAC-2XMO	Positive, 2X Cutoff with MDMA and Oxycodone	5 mL
DOAC-19572753	GC/MS-H Cutoff+25%	50 mL	DOAC-2XMOB	Stat-Skreen, 2X Cutoff with MDMA, OXY & BUP	5 mL
Confirm 4 HIGH OPIATE			DOAC-F	Positive, Fentanyl (norfentanyl)	5 mL
DOAC-19582041	Confirm 4, Cutoff -60%	50 mL	NEGATIVE		
DOAC-19582026	Confirm 4, Cutoff -25%	50 mL			
DOAC-19582002	Confirm 4, Cutoff	50 mL	DOAC-19227000	Negative Control Urine	50 mL
DOAC-19582126	Confirm 4, Cutoff +25%	50mL	DOAC-000-1	Negative Control Urine	5 mL
Confirm 5 HIGH OPIATE			GC/MS MIX		
DOAC-18008614	Confirm 5, Cutoff -60%	50mL	DOAC-18002988	GC/MS Mix	20 mL
DOAC-18008615	Confirm 5, Cutoff +25%	50mL			

Nova Century Scientific, Inc.
A Division of Trinity Biotech Plc.
5022 South Service Road
Burlington, Ontario L7L5Y7, Canada

Phone: 1 800 615 5072
 Fax: 1 800 639 9006

Email: nova.info@trinitybiotech.com
 Website: www.trinitybiotech.com

For additional information on our other products please contact us or refer to our website.

SYMBOL LEGEND	
	Consult Instructions for Use
	Temperature Limits
	In Vitro Diagnostic Medical Device
	Batch Code
	Product Catalog Number
	Manufacturers Identification
	Use by Date
	Caution, Consult Accompanying Documents

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